


For Any Technical Questions, Please Contact at info@nodesus.com

Product Name: TRANSFORM Microwave Digestion System

Part Number:

Product Description:

TRANSFORM Microwave Digestion System(TRANSFORM 680/TRANSFORM 800):



TRANSFORM 680 microwave digestion systems



TRANSFORM 800 microwave digestion systems

Aurora's TRANSFORM Series of automatic top-loading, closed-vessel microwave digestions systems are equipped with high-pressure capabilities to enhance digestion quality while decreasing digestion time. These systems can simultaneously run up to 10 high-pressure vessels, while providing a fast, automated method to digest even the most difficult samples.

The TRANSFORM Series provides unparalleled safety, reproducible measurements, and ease of operation, allowing you to streamline your day-to-day operations. With two possible capacities (six or ten sample), up to 1200 W of power output, and continuous pressure and/or temperature monitoring, our TRANSFORM MW units can meet your digestion needs. These systems prepare samples for a range of downstream analyses, including AAS, AFS, Hg analyzers, ICP-OES, and ICP-MS.

Performance Features



TRANSFORM 680 microwave digestion systems

- √ Simultaneous loading and unloading of six (6) high pressure digestion tanks (up to 800 PSI)
- √ Direct measurement of sample temperature with a precision of +1°C
- √ Suitable for food, metallurgy, and electronics digestion



TRANSFORM 800 microwave digestion systems

- √ Simultaneous loading and unloading of ten (10) high pressure digestion tanks (up to 800 PSI)
- √ Real-time monitoring of temperature and pressure during digestion
- √ Top-loading auto-rotating door provides heavy-duty protection
- √ More suitable for the treatment of petroleum, plastics, and other substances

For Any Technical Questions, Please Contact at info@nodesus.com

Analytical Performance

- √ Shortens digestion time (when compared with traditional digestion methods)
- √ Reduces loss of volatile analytes and components
- √ Eliminates environmental sample contamination
- √ Reduces solvent consumption
- √ Allows for unattended operation
- √ Online EPA Method available
- √ In-situ air cooling technology: fast, safe, and efficient

Standards and Certification

- √ Compliance with RoHS and WEEE standards
- √ Passed ISO 9001 and CE certification

Instrument Compliance

The TRANSFORM series conforms to various EPA methods and ASTM related standards. For user's convenience, EPA method has been pre-installed in the software. Any of these methods can be utilized to accommodate different types of sample preparation.

- √ EPA 3015
Microwave-assisted Digestion of Aqueous Solution Samples
- √ EPA 3051
Microwave-assisted Digestion of Sediments, Mud, Soil, and Petroleum
- √ EPA 3052
Microwave-assisted Digestion of Organic Samples and Silica in Soil


nodesus
Engineering Your Needs

Features and Performance



TRANSFORM 800 Top View




TRANSFORM 800 Digestion Tank



Main Control Tank Connections

Software Interface



- ✓ Safe external control via PC and software, allows remote control of magnetron, cooling fan, sample rotating disc, etc.
- ✓ Real-time display of temperature and pressure curves during digestion
- ✓ Real-time display of microwave emission power and temperature change of magnetron during digestion
- ✓ Easy to change sample digestion procedures, such as adding/deleting digestion procedures
- ✓ Pre-installation of EPA methods and parameters makes the digestion of various samples more convenient

Heavy Top Loading Furnace Body Safety Features

- Top-loading automatic door provides heavy-duty blast protection
- Double chamber, explosion-proof safety door, with inner chamber door spring device that can effectively absorb explosion shockwave and reduce blast pressure
- Automatic (computer-controlled) safety interlocking device secures blast door, whilst providing instrument insulation that ensures microwave leakage is $< 5\text{mW/cm}^2$

Automatic Door Safety Features

- Opening and closing of the door is computer controlled, reducing staff exposure to explosion risks
- Software continuously monitors the opening/closing process, ensuring that the door is closed before instrument sequence can begin

Safety Tank Features

- Double Layer Tank: High strength PEEK plastic-coated external tank and thickened TFM material
- Internal Tank: High temperature and pressure-tested, corrosion resistant
- Explosion-prevention film: film ruptures when pressure exceeds safety value, minimizing over-pressure and explosion risk



Outer tank



Inner tank

Preventive measures for acidic vapors

- Three-layer Teflon Antiseptic Surface Treatment
- Acidic Vapor Tank: collects acidic vapors released by the rupture of rupture-disks on the digestion tanks, eliminating acidic vapor diffusion, thus preventing furnace body corrosion and prolonging instrument service life
- Exhaust and Cooling: quickly removes acidic vapors from the chamber, followed by rapid air-cooling of samples

Precise control of the digestion process

- Platinum-sensors, monitor temperature and pressure; highly sensitive, guarantees accurate monitoring and control of temperature and pressure
- Real-time visual monitoring of temperature and pressure to ensure actual temperature consistent with preset temperature curve
- The sample disc rotates continuously in one direction, ensuring a more uniform distribution of heat
- Microwave non-pulse continuous emission, allowing microwave intensity adjustments according to the degree of heating that is required

Pressure Sensor



Temperature Sensor

For Any Technical Questions, Please Contact at info@nodesus.com

Product Specification:

| | TRANSFORM 800 | TRANSFORM 680 |
|------------------------------------|--|---|
| External Dimension(mm) | (L)550×(W)490×(H)490 | (L)480×(W)430×(H)430 |
| Power Source | 220V, 50Hz, 12A | 220V, 50Hz, 12A |
| Weight (kg) | 50 | 36 |
| Microwave Power (W) | 1200 | 1200 |
| Cavity Volume (L) | 34 | 18 |
| Explosion-proof safety door | Included | Included |
| Sample Loading method | Top Loading, Electronic Control Safety Door | Top Loading, Electronic Control Safety Door |
| Temperature Control System | Platinum-coated, corrosion-resistance probe in reference tank, able to withstand wider range of temperatures, providing high accuracy. | |
| Pressure Control System | Reference tank piezoelectric crystal pressure sensor | N/A |
| Sample positions | 10 | |
| Cooling Mode | Air Cooling | Air Cooling |
| Sample Reaction Tank: Outer Tank | PEEK | PEEK |
| Sample Reaction Tank: Inner Tank | TFM | TFM |
| Volume of each digestion tank (mL) | 50 | 80 |
| Maximum Pressure (Psi) | 2000 | 1600 |
| Critical Temperature(°C) | 300 | 300 |
| Maximum Work Pressure (psi) | 800 | 400 |
| Maximum operating temperature(°C) | 250 | 250 |