


Deck Layouts



VERSA Sample (Environmental) Processing Deck
Aspirates environmental samples from tall tubes and dispenses into specially designed recipient racks or carousels. Capable of user-defined dilutions, repetitive and complex environmental liquid handling processes. Sonication washing of the pin after every transfer minimizes contamination of samples.



VERSA Mini Vial Filler Deck
Liquid dispensing system for vials and tubes (48 and 96 tube formats).



VERSA Cap Piercing Deck
4 - channel cap piercing head automatically pierces sample container, aspirates sample, and dispenses into target plate.



VERSA Mini PCR Deck-
Shown with Roche LightCycler & GeneDisc AdaptorSets
VERSA automated liquid handling workstations compatible with Roche's light cycler carousel adaptor and Corbett's GeneDisc Adaptor.



VERSA Mini PCR Deck-
Shown with RG 6000 & GeneDisc AdaptorSets
VERSA automated liquid handling workstations compatible with Qiagen's RG-6000 and GeneDisc adaptor



VERSA Mini PCR Deck-
Shown with RG 6000 & Roche-Cobas AdaptorSets
VERSA automated liquid handling workstations compatible with Qiagen's RG-6000 and Roche-Cobas Ring Adaptor.